

HBG



HBG SAHLI'S HAEMOMETER

The HBG™ haemometer as a complete set consists of:

- Polystyrene support with 2 colored rods and opal glass plate.
- Comparator tube.
- Hemoglobin pipette 20 μ l acc. to Sahli.
- Silicone tubing of approx. 16 cm length, white mouthpiece.
- Dropping pipette with rubber teat.
- Stirring rod, Acid vial.
- Cleaning brush, Directions for use (Manual).

Pipetting by mouth is dangerous. Please ensure safety while using it. Read the instruction manual carefully.