

GLOLAB

I N S T R U M E N T S

Tablet Disintegration Apparatus GDA-101



FEATURES

1. Complies with USP,IP Specifications
2. Extremely useful for Pharmaceuticals labs
3. Highly Accurate and easy to Operate
4. Baskets 'UP' & 'DOWN' Movement as per IP/USP Specifications
5. 20 x 4 line Alphanumeric LCD with Backlit
6. Printer Attachment Facility
7. Audio Visual Indicators for Systems Status.

TECHNICAL SPECIFICATION

Display	: 20 x 4 Line back lighted LCD Display
Keyboard	: Alphanumeric Function Keys
Dip Speed	: 30, ± 1 DPM
Stroke Length	: 55, ± 1 mm
Beaker Volume	: 1000 ml x 2 Beakers
Basket	: 2 x 6 Test Positions (10MESH) USP type A
Fluted Disk	: 6 Numbers in Each Basket USP type A
Water Bath Temp.Range	: Ambient +5 to 40°C
Water Bath Temp.Accuracy	: ± 0.3 °C
Water Bath Circulation	: Magnetic Pump
Test Mode	: Programmable Time up to 99:59(HH:MM) Display Format 99:59:59
Output Printer	: Serial Port
RS-232C	: PC Connectivity for Data Downloading
Dimension (Approx.) [LXBXH]	: 420 X 380 X 520 mm
Shipping Weight (Approx.)	: 18 kg
Electrical Power	: 230V \pm 23V, 50/60 Hz or 115V \pm 15V, 50/60 Hz