Rico Scientific Industries

- **Automatic Calculation & Display of Calorific Value.**
- **Peak Temp. Hold with Alarm.**
- **❖** With Inbuilt Printer. Calorific Value, Date & Time is printed on the Report Generated.
- **❖** Readability of 0.01 °C in Ambient & Rising Temp.
- **Accuracy of ± 0.01°C in Rising Temp.**
- **Accuracy of Calorific Value ± 50 calories.**
- Menu Key is provided to select type of Test i.e Calibration / Water Equivalent, GCV Sample Test.
- ❖ Numeric Digits 0-9 are provided on the front panel for feeding Values such as Weight of the tablet, Water Equivalent values etc.
- ***** Waterproof Long Life Brass Sensor.

SUPPLIED COMPLETE WITH 26 ACCESSORIES LIKE PUFF-INSULATED JACKET, SAFETY DEVICE, CERTIFIED S.S. BOMB (BY GOVT. OF INDIA), TABLET MAKING MACHINE, OXYGEN REGULATOR WITH PRESSURE GAUGE, SAETY DEVICE & FLEXIBLE TUBING, SPARES ETC.

PRE INSTALLATION REQUIREMENT:- FRESH OXYGEN CYLINDER (IND<u>USTRIAL GRADE)</u>, DIGITAL WEIGHING BALANCE ATLEAST (0.001G), 5-10 LITERS DM / DISTILLED WATER.

RICO Microprocessor Based Automatic Bomb Calorimeter Model RSBT - 7



PRINT FORMAT:-

BOMB CALORIMETER

WATER EQUIVALENT TEST

MAX.TEMPERATURE 02.67 C

BENZOIC ACID CV 6300

MASS OF TABLET 1.000 GM

WATER EQUIVALENT 2371

18:15 19/09/13

BOMB CALORIMETER

SAMPLE TEST

MAX.TEMPERATURE 00.05 C

WATEREQUIVALENT 2500

MASS OF TABLET 0.995 GM

CV 95.32 CAL/GM

14:34 09/09/13