

HALOGEN DRYING METHOD

MOISTURE ANALYZER LH

76

NEW



TECHNICAL SPECIFICATIONS

✓ Weighing Range	50g
✓ Readability	1mg
✓ Redability (3g samples)	0,20%
✓ Minimum Sample amount	0,5g
✓ Advised Sample amount	3-10g
✓ Heat up time	1-99 minutes, 1 minute interval
✓ Temperature Program	Standard
✓ Terminal Control	Timing, Automatic
✓ Heating Temperature Range	50-180°C
✓ Showing Content	Moisture %, Solid %, Weight, Time, Data etc.
✓ Heating Source	Halogen Lamp
✓ Pan Size	100mm
✓ Dimensions	265 x 160 x 150 mm
✓ Net Weight	5 Kg
✓ Package Size	530 x 380 x 340 mm
✓ Gross Weight	6 Kg



CE FDA

Laboholic

31